
USER MANUAL

HELIX Mobility Scooter

Model: MJMA01

By Travel Buggy

★CAUTION: Read all precautions and instructions in this manual before using this equipment.

Improper use or maintenance can void the warranty.
Keep this manual for future reference.

NOTICE

Thank you for choosing HELIX Mobility Scooter by Travel Buggy, a convenient to use and foldable 4 wheel steering power scooter.

It is important to read through this manual before operating your HELIX Mobility Scooter to become familiar with its features as well as its limitations and safety information

Operating and maintaining your HELIX Mobility Scooter is designed to be straightforward and hassle-free. However, if you have any queries or need assistance, feel free to reach out to us at Travel Buggy.

OUR CUSTOMER HOTLINE

1.888.56.BUGGY (Canada)

1.855.56.BUGGY (USA)

EMAIL : contact@travelbuggy.com

DISCLAIMER:

The notes, warnings, and cautions contained within the manual and marked by the triangular Caution Symbol should be given special care. Users should also pay special attention to information marked in this manual beginning with NOTICE.

Travel Buggy, its affiliates and partners cannot be held responsible for personal injury or property damage resulting from the unsafe or improper operation or maintenance of the HELIX Mobility Scooter.

As part of our ongoing product improvement initiative, Travel Buggy reserves the right to change specifications and design without notice.

As a result, there may be minor differences between your HELIX Mobility Scooter and the photos, illustrations and instructions in this manual.

SAFETY INFORMATION

- DO make certain that the seat is in a locked position when you are operating your HELIX
- DO reduce your HELIX' s speed when you are driving around a corner
- DO proceed with caution as you approach the downgrade of a ramp or other incline and when driving near a raised surface or unprotected ledges and drop offs (curbs, porches, stairs, etc.)
- DO avoid exposure to rain, snow, ice, salt, or standing water whenever possible
- DO maintain and store in a clean and dry condition
- DO keep your hands away from the tires when driving. Be aware that loose fitting clothing can become caught in drive tires
- DO keep tools and other metal objects away from battery terminals. Contact with tools can cause electrical shock
- DO follow the recommended weight capacity (264 lbs / 120 kg)
- DO follow the recommended incline slope rating (12 degrees)
- DO check the battery charge to ensure there is enough charge to get to your destination before using.
- DO ensure your HELIX' s power is OFF when entering and exiting the power mobility scooter
- DO NOT drive your HELIX across the side of a hill or diagonally up or down a hill, and do not stop (if possible) while driving up or down an incline
- DO NOT operate your new HELIX for the first time without completely reading and understanding the owner' s manual
- DO NOT carry additional passengers on your HELIX
- DO NOT allow unsupervised children to play near the power scooter
- DO NOT operate your HELIX while under the influence of alcohol, as this may impair your ability to drive safely
- ALL servicing must be performed by qualified personnel only. Unauthorized modifications to power scooter, its power system or its control system may compromise the performance and safety of the system. UNDER NO CIRCUMSTANCES should you attempt any repair or service that is not covered in this manual.
- ALWAYS ensure all safety warning labels are in place and legible.

BATTERY SAFETY INFORMATION

1. Lithium-ion battery packs may get hot, explode or ignite and cause serious injury if exposed to abuse conditions. Be sure to follow the safety warnings listed below:

- DO NOT place the battery near fire or heat the battery
- DO NOT install the battery backwards so the polarity is reversed
- DO NOT connect the positive and negative terminal of the battery to each other with any metal object (such as a wire)
- DO NOT carry or store the battery together with necklaces, hairpins, or other metal objects
- DO NOT pierce the battery with nails strike the battery with a hammer, step on the battery or otherwise subject it to strong impacts or shocks
- DO NOT solder directly onto the battery
- DO NOT expose battery to water or salt water, or allow the battery to get wet

2. DO NOT disassemble or modify the battery. The battery contains safety and protection devices, which if damaged, may cause the battery to generate heat, explode or ignite

3. DO NOT place the battery in or near fire, on stoves or other high temperature locations. Do not place the battery in direct sunlight, or use or store battery inside cars in hot weather. Doing so may cause the battery to generate heat, explode or ignite. Using the battery in this manner may also result in a loss of performance and a shortened life expectancy.

4. If the device is to be used by children, the caregiver should explain the contents of this document to the children and provide adequate supervision to ensure the device is being used appropriately

5. When the battery is worn out, insulate the terminals with adhesive tape or similar materials before disposal

6. Immediately discontinue use of the battery if, while using, charging or storing the battery, the battery emits an unusual smell, feels hot, changes color, or shape, or appears abnormal in any way. Contact the Fire Department if any of these problems are observed

7. DO NOT place battery in microwave oven, high pressure containers or on induction cookware

WARNING:

The operator should read all of the warning signs and instructions here and on the power scooter before use.

You should not operate your power scooter on public streets and roadways. Be aware that it may be difficult for traffic to see you when you are seated on your power scooter. Obey all local pedestrian traffic rules. Wait until your path is clear of traffic, and then proceed with extreme caution.

BATTERY SAFETY INFORMATION

CHARGING PRECAUTIONS

1. BE SURE to follow the warnings listed below while charging the battery. Failure to do so may cause the battery to become hot, explode or ignite and cause serious injury

- When charging the battery, ONLY USE THE SPECIFIED CHARGER PROVIDED WITH THE PRODUCT

- DO NOT attach the battery to a power supply plug or directly to a car's cigarette lighter

- DO NOT place the battery in or near fire, or in direct sunlight. Heating the battery can damage the safety circuitry, which can cause additional heating, rupture or ignition of the battery

2. DO NOT continue charging the battery if it does not recharge in the timelines specified in the operating manual

3. THE TEMPERATURE range over which a battery can be charged is 5°C to 40°C (41°F to 104°F). Charging the battery at temperatures outside of that range may cause severe damage to the battery or reduce battery life expectancy

- Always protect battery from freezing temperatures and never charge a frozen battery. This will damage the battery and may cause personal injury. Attempting to charge a battery in freezing conditions does not prevent a battery from freezing

WHEN DISCHARGING THE BATTERY

1. DO NOT use the battery with any device other than the HELIX. Utilizing the battery in devices other than the HELIX can lead to damage or a reduction in its lifespan.



Additionally, abnormal current flow caused by other devices may result in the battery overheating, potentially leading to explosion, ignition, and serious injury

2. THE TEMPERATURE range over which the battery can be discharged is -20°C to 60°C (-4°F to 140°F). Use of the battery outside this temperature range may damage the performance of the battery or may reduce its life expectancy

TABLE OF CONTENTS

Main Components	6
Specification	7
Operation Instructions	8
★ Unboxing instruction	8
★ Handrail installation (Optional	10
★ Power switch	11
★ Speed adjustment	11
★ Driving Control	11
Lights and horns	12
Brake lever	12
Folding method	12
Battery and charging	13
★ Battery considerations	13
★ Charging method	13
★ Charger	14
Motor Controller and Connections	14
Maintenance	15
Daily inspection	15
Routine Maintenance	15
Troubleshooting	17
Limited Warranty	18
Contact & Support	19

MAIN COMPONENT

- | | |
|-----------------------------------|-----------------------------|
| ① Power switch & speed adjustment | ⑩ Detachable seat pad |
| ② Head light button | ⑪ Battery |
| ③ Throttle & brake | ⑫ Seat lock operating lever |
| ④ Height adjustment latch | ⑬ Anti tip wheel |
| ⑤ Front Angle adjustment lever | ⑭ Taillight |
| ⑥ Gas shock | ⑮ Motor lock |
| ⑦ Horn button | ⑯ Front trim |
| ⑧ Footrest | ⑰ Head light |
| ⑨ Backrest pad | ⑱ Megaphone |

HELIX Mobility Scooter

SPECIFICATIONS

Package size (L*W*H)	920*640*380mm
Unfold size (L*W*H)	1020*610*350mm
Fold Size (L*W*H)	890*610*350mm
Seat cushion size (L*W*H)	400*350*70mm
Gross weight	29.8KG(excluding armrests)
Gross weight(excl battery)	27.4kg
Seat height from the ground	465mm
Handlebar height from the ground	800-850mm adjustable
Maximum load	120kg
Motor	250W brushless
Battery	Lithium battery 24V 5.8Ah * 2
Controller	Intelligent controller
Headlights	LED
Taillight	LED
Screen	LCD instrument
Rear wheel size	8 "PU solid tires
Front wheel size	8 "PU solid tires
Speed	0 ~ 6 km/h adjustable
Climbing Angle	≤ 12°
Climbing distance	≤ 200m
Steering	Rear wheel follow-up steering
Minimum turning radius	97cm
Obstacle crossing height	≤ 5cm
Travel distance	20km
Brake	Electromagnetic brake
Charger	Special DC24V/2A AC110-250V
Charging time	4-6h

7. TRAVEL BUGGY

OPERATION INSTRUCTIONS

★ UNBOXING

1. Find the upward symbol “ ↑ ↑ ” on the outer packaging box and place the carton in the direction of the arrow. First, open the carton from the top and observe the position of the scooter's wheels. Remove the outer carton and take out the scooter accessories, pedal accessories and excess sponge MATS. Observe whether there are bumps in the scooter and whether the accessories are complete.



Please ask for assistance if necessary, as wrong operation may lead to damaging the mobility scooter.

2. First pull down the knob at both ends of the pedal cross bar and turn it to the bump. Take out the pedal fittings and align them with the circular holes at both ends of the bar.

3. Insert the U-shaped pipe on the back of the foot pedal into the circular holes at both ends of the cross bar and restore the knob so that the top of the knob is stuck on the foot pedal. At the same time, the two ends of the "click" sound is heard, and the foot pedal is installed.


HELIX Mobility Scooter


4. Press the red button with one hand and press the handle with the other hand, slowly lifting it to the appropriate angle.

Gas shock is under pressure, please make sure your both hands to operate it at the same time, when extending the gas shock.




5. Then, lift the backrest. Pull the release lever upward with one hand while slowly holding the seat with the other, when you hear a click, release the lever and let it fall. At this point, the seat should remain level, and then place the seat cushion. Otherwise, repeat the above steps until level.

6. Open the handlebar height adjustment buckle, you can adjust the handlebar height up and down, and then lock the buckle. Finally, place the battery lock upward, place the battery in the battery carrier until it locks.

HELIX Mobility Scooter

★ Handrail installation (Optional)

1. First open the knobs at both ends counter clockwise, and then insert the armrest into the square tube. Adjust the appropriate armrest position and tighten the knob clockwise to fix the armrest position.
2. Adjust the backrest by squeezing the backrest folding switch inwards.

★ Check the mobility scooter, footrest and battery .etc. for damage. (if there is any damage, please contact the dealer for replacement)

★ Power switch

- 1 This scooter can be equipped with two battery packs, the battery has a separate switch located at the front of each battery box, as shown in the figure, "○" indicates that the power off, "—" indicates that the power on.
- 2 Make sure that the battery power has turned on, rotate the round knob on the control panel clockwise until you hear a "click", then the screen will lite up, now the mobility scooter has been turned on and you can drive
- 3 Depending on your needs, you can use one battery separately or use two batteries at the same time.

Please do not turn off the power when the mobility scooter is driving, otherwise the electromagnetic brake will engage and the scooter will stop suddenly.

★ Speed adjustment

The maximum forward speed can be adjusted within the range of 0.2Km/h ~ 6Km/h, which can be set. Rotate round knob counterclockwise for a slow speed setting, rotate the round knob clockwise for a fast speed setting.

★ Driving Control

- 1 Pressing the right paddle on the throttle will make the scooter go forward. Pressing the left throttle paddle will make the scooter go in reverse.

- 2 After the throttle is released, the scooter will automatically stop in the center position. The electric brake will engage and stop the scooter.

HELIX Mobility Scooter

★ Lights and horns

Pressing the yellow button will sound an alarm; Press the red button to turn on the headlight and press again to turn it off.

★ Brake lever

A red lever under the taillight can be switched between drive and a manual mode.
"Drive": Push the brake lever down (Locked) for mobility mode.
"Manual mode": Pull the brake lever up (unlocked) for manual mode.

- Please turn off your scooter when you are done driving.
- Do not press the throttle with both hands at the same time.
- Do not turn off the power switch while driving, otherwise, this will lead to an emergency stop and can cause you to fall off the scooter.
- The braking distance will change according to the driving road conditions, so if you want to stop, please brake the mobility scooter in advance.

★ Folding method

1
Switch to manual drive mode.

2
Remove the seat cushion and backrest cushion and place them in the appropriate position.

3
Lift the folding lever with one hand and support the seat with the other hand, slowly lowering the seat.

4
Press the red button with one hand and the handlebar with the other hand to fully fold it.


5
Finally, the backrest is folded.

BATTERY AND CHARGING

★ Battery informations

This product uses two sets of 24V 5.8AH high-performance lithium-ion batteries, which must be fully charged for the first use.

- ① Power
- ② Charging port
- ③ USB
- ④ Battery release switch
- ⑤ Battery status indicator

- Do not disassemble the lithium battery privately, otherwise it may cause accidents and body injury.
- Battery operating temperature: $5 \pm 2^{\circ}\text{C} \sim 40 \pm 2^{\circ}\text{C}$.
- Don't drive scooter on low battery capacity for a long time, otherwise the battery life will be shortened.
- Please charge the battery every two months, charge it about 80% of battery capacity. If the battery loses power for a long time, the battery will be damaged, and will void your warranty.

★ Charging method

Switch to manual mode.

Plug the charger into the charging port of the battery under the seat.

Insert the charger power plug into the home power supply.

When the charger is connected to the wall power supply, the indicator lights up red, and the charging time according to the battery state and temperature which is generally 2-6 hours.

After charging, the charger automatically stops working and the indicator changes to green, please pull out the input plug of the charger, and then pull out the plug of the battery side.



★ Charger

- 1.This product uses AC100-240V / DC24 / 2A lithium battery special charger, do not mix other chargers as this will void your warranty.
- 2.Do not disassemble or modify the charger.
3. Do not cover and block the charger to avoid affecting heat dissipation
4. Keep the charging environment ventilated.

- When operating the scooter, both batteries will discharge simultaneously .
- When charging using the in-seat port, both batteries will charge simultaneously.
- Do not charge using the in-seat port when one battery is discharged, and one is full, this will shorten the battery life.
- The scooter can be operated using only one battery, either the front or rear battery slot.

★ Motor, Controller and its connections

1

2

3

4

1. Controller connection to motor
2. Controller connection to battery
3. Controller connection to rear LED light
4. Controller connection to operation panel and front LED light
5. Controller unit

★ Daily Inspection

Inspection item	Scope of examination
Handle	Is there a loose, whether it can be normal around? Is the expansion normal?
Seat and backrest	Is there any loose, inclined, and whether it can fold normally?
Speed adjustment knob	Does it rotate normally?
Forward and backward control lever	Can you adjust the speed normally?
	Can you push it down and return to the middle and stop?
1 Brake lever	Normal separation?
3 Battery status	Does the battery meter match the control panel? Is the battery charge sufficient?

★ Routine Maintenance

Check regularly insulation of all wires and tightening of wire end regularly.

Avoid moisture of all electronic control module, otherwise it must be completely dry before use.

Please use a soft duster cloth to wipe the dirt of the mobility scooter body and do not clean the scooter body directly with water.

TROUBLESHOOTING

The controller contains a program to identify problems and troubleshoot. When the diagnostic program detects problems and errors, the LED on the control panel shines, the fault alarm buzzer will alarm, and the fault display LED shines continuously in the form of a 1-bit fault indicator code until troubleshooting.

If the problem persists, please get in touch with us or the dealer.

Error code	Fault description	Failure category	Solutions
E0.1	Over current fault	Drive failure	a)Check whether the control panel line sequence is correct b)Is the fault of the replacement instrument normal c)Is it normal to replace the controller d)Replace the instrument and controller at the same time
E0.2	Under voltage fault	Drive failure	a)Check battery power and charge b)Trouble shooting replacement of drive
E0.3	Over voltage fault	Drive failure	a)Check battery power b)Troubleshooting replacement drive
E0.4	Stalled rotor fault	Drive failure	a)Check whether the motor is stuck b)Troubleshooting replacement of drive c)Troubleshooting replacement of motor
E0.5	The above mos fault	Drive failure	a)Troubleshooting replacement of driver b)Check the motor replacement
E0.6	The following mos fault	Drive failure	a)Troubleshooting replacement of driver b)Check the motor replacement
E0.7	Hall failure	Drive failure	a)Check whether the motor wiring harness is damaged, disconnected or short-circuited b)Troubleshooting replacement of drive c)Check the replacement motor

TROUBLESHOOTING

Error code	Fault description	Failure category	Solutions
E0.8	mos over temperature fault	Drive failure	a)Cool down the drive b)Troubleshooting replacement of drive c)Check the replacement motor
E1.2	communication fail	Instrument/ drive failure	a) Check whether the instrument line sequence is correct b) Is the fault of the replacement instrument normal? c) Is it normal to replace the controller? d) Replace the instrument and controller at the same time
E1.7	Brake failure	Drive failure	a)Check whether the motor brake is closed b)Check whether the motor wiring hamess is damaged or short-circuited c)Troubleshooting replacement of drive d)Check the motor replacement
E2.1	Forward and backward potentiometer falls off	Instrument failure	a) Check whether the handle is well connected b) Check whether the handlebar line sequence is correct
E2.2	The forward and backward potentiometers are not in the neutral position when turned on	Instrument failure	a) Check whether the handlebar is in the neutral position b) Check whether the handlebar line sequence is correct

HELIX Mobility Scooter

LIMITED WARRANTY

The purchaser of this product is entitled to a limited warranty as offered by Travel Buggy and its affiliates on the following components and timeframes from the date of purchase:

Warranty Component	Warranty period
Main frame of mobility scooter	60 Months
Motor	12 Months
Battery	12 Months
Controller	12 Months
Wear Parts & Charger	3 Months

★ This warranty is valid for the replacement of dysfunctional parts only. Any parts under warranty will be replaced and shipped for free. Any service and labor fees, if applicable, to replace parts under warranty must be paid by the user.

Due to its straightforward design, most parts can be easily exchanged by the end user without a professional service tech required. However, it is always recommended you seek professional help for maintenance and service, to make sure the work is done properly. The warranty only covers NON-DELIBERATE Damage.

★ The warranty does not cover:

- Products damaged by user negligence.
- Products damaged accidentally.
- Products damaged intentionally.
- Products that have been subjected to negligence.
- Products that have been subjected to abuse.
- Products that have been improperly stored.
- Products that have been improperly handled.
- Products that have been improperly operated.
- Products that have experienced general misuse.
- Products that have been modified in an unauthorized, unapproved way.

★ SERIAL NUMBER

The serial number on the HELIX is located on the front tiller underneath the rubber cover . It is important to mention the serial number in any warranty related communications

Warranty is non-transferable and only valid for the original scooter purchaser. Travel Buggy will, at its sole discretion, replace any warrantable item that fails within the warranty period.

CONTACT & SUPPORT

Operating and maintaining your HELIX is designed to be straightforward and hassle-free. However, if you have any queries or need assistance, feel free to reach out to us at Travel Buggy.

OUR CUSTOMER HOTLINE

★ 1.888.56.BUGGY (Canada)

★ 1.855.56.BUGGY (USA)

EMAIL

★ contact@travelbuggy.com


Additional images and videos are available on our website.

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