



BOLERO

Sturdy & Rugged

Hoveround's Bolero is the ideal heavy duty mobility scooter for outdoor use. Its 350-pound weight capacity and deluxe captain's seat makes it one of the most sturdy, durable and stable 4-wheel electric scooters of its kind.

FEATURES & BENEFITS

- Rugged, four-wheel design makes for a sturdy and stable ride.
- Dual rear-view mirrors, headlights and rear lights for added safety.
- Adjustable, flip-up armrests for comfort and safe transfer to scooter.
- Deluxe 360-degree swivel captain's seat and headrest provide exceptional riding comfort.
- Attached basket for transporting items while on-the-go.
- Lightweight, portable charger lets you charge your batteries anywhere.

SPECIFICATIONS

Max weight capacity.....	350 lbs
Max speed ¹	6.3 mph
Max range per charge ¹	20 miles
Turning radius.....	43.3"
Ground clearance	3"
Length.....	51.5"
Width	23.6"
Drive tires	(2) 10" x 3.3"
Front tires.....	(2) 10" x 3.3"
Total weight with batteries	195 lbs

¹ Speed and range will vary with respect to terrain and weight of rider. The Manufacturer reserves the right to alter specifications without notice.

HOVEROUND[®]

Personal Mobility Vehicles