



HOVEROUND MPV5

Most maneuverable

The MPV5[®] is not just one of the most maneuverable power chairs available today, it is also Hoveround's most popular model. The unique, Round for a Reason[®] design enables easy navigation around tight corners and through narrow spaces without damaging walls or furniture. The MPV5 can safely function indoors as well as outdoors, offering power wheelchair users greater flexibility in their freedom and mobility.

FEATURES & BENEFITS

- Flip-up footplate, round frame and compact design for easy maneuverability.
- Custom-contoured swivel seat for utmost comfort and ease of transfer.
- Sturdy cup holder and side tote bag for increased independence and convenience.
- Adjustable armrests for comfort and maximized support.
- Lightweight, portable charger for recharging power chair batteries anywhere.

SPECIFICATIONS

Max weight capacity.....	300 lbs
Max speed ¹	5 mph
Max range per charge ¹	15-20 miles
Turning radius.....	22.7"
Ground clearance	2.5"
Maximum step climb	2"

Length.....	38"
Width	24"

Drive tires	(2) 9" x 3.5"
Caster tires	(2) 8" x 2"

Total weight with batteries	195 lbs
Base only.....	91 lbs
Base with batteries.....	143 lbs

¹ Speed and range will vary with respect to terrain and weight of rider. The Manufacturer reserves the right to alter specifications without notice.

