



TEKNIQUE FWD

Front-wheel drive

The Hoveround Teknique® FWD power chair is designed to easily maneuver both in the home and outdoors. Everyone will appreciate the comfortable design and powerful front-wheel drive maneuverability. With automatic braking to ensure safe vehicle operation and seven base colors to choose from to make it yours, the Teknique FWD makes a great choice for improved mobility.

FEATURES & BENEFITS

- The better choice for individuals over 6 feet tall.
- Flip-up footplate, round frame and compact design for easy maneuverability.
- Custom base plate available in black, red, teal, purple, yellow, blue or white.
- Higher ground clearance and powerful front-wheel drive performance.
- Sturdy cup holder and side tote bag for increased independence and convenience.
- Adjustable armrests for comfort and maximized support.
- Lightweight, portable charger for recharging power chair batteries anywhere.

SPECIFICATIONS

| | |
|---|----------------|
| Max weight capacity..... | 300 lbs |
| Max speed ¹ | 4.5 mph |
| Max range per charge ¹ | 15-20 miles |
| Turning radius..... | 24" |
| Ground clearance | 3.6" |
| Maximum step climb | 2" |
| Length..... | 41.9" |
| Width | 25" |
| Drive tires | (2) 10" x 3.5" |
| Caster tires | (2) 8" x 2" |
| Total weight with batteries | 179 lbs |
| Base only..... | 90 lbs |
| Base with batteries..... | 139 lbs |

¹ Speed and range will vary with respect to terrain and weight of rider. The Manufacturer reserves the right to alter specifications without notice.

