



TEKNIQUE GT

For active users

The Teknique® GT power chair combines the sturdy design and custom-contoured seating, found throughout the Hoveround Teknique line, with additional power for the active user. Strong enough to safely maneuver outdoors, this rear-wheel drive power chair features larger drive tires that give more active users superior riding performance with the expanded mobility they seek.

FEATURES & BENEFITS

- Fixed-angle, swing-away footrests with heel loops for secure leg support.
- Large drive tires increase ground clearance for maneuvering over inclines and small obstacles
- Sturdy cup holder and side tote bag for increased independence and convenience.
- Adjustable armrests for comfort and proper upper body support.
- Custom-contoured seat for maximum comfort.
- Dual motors with powerful rear-wheel drive for extra riding stability.
- Lightweight, portable charger for recharging power chair batteries anywhere.

SPECIFICATIONS

Max weight capacity.....	300 lbs
Max speed ¹	6.5 mph
Max range per charge ¹	15-20 miles
Turning radius.....	37"
Ground clearance	4.7"
Maximum step climb	3"
Length.....	42.5"
Width	25.3"
Drive tires	(2) 13" x 4"
Caster tires	(2) 9" x 2.8"
Total weight with batteries	242 lbs
Base only.....	102 lbs
Base with batteries.....	181 lbs

¹ Speed and range will vary with respect to terrain and weight of rider. The Manufacturer reserves the right to alter specifications without notice.

