



TEKNIQUE XHD

Bariatric support

The highly maneuverable Technique® XHD is the definitive bariatric power chair for anyone in need of a more robust mobility solution. The XHD is designed with a 450-pound weight capacity and runs on heavy-duty 70 amp batteries for extended use. Combined with wider seating for comfort and tandem rear tires for increased traction, the Technique XHD stands out as the bariatric power chair of choice in its class.

FEATURES & BENEFITS

- Durable, heavy-duty construction for bariatric users from 300 - 450 lbs.
- Adjustable height and back positioning for extra comfort.
- Custom-contoured wide seat (22" – 24") for maximum riding support and comfort.
- Sturdy cup holder and side tote bag for increased independence and convenience.
- Adjustable armrests for comfort and proper upper body support.
- Oversized footplate.
- Lightweight, portable charger for recharging power chair batteries anywhere.

SPECIFICATIONS

Max weight capacity.....	450 lbs
Max speed ¹	5 mph
Max range per charge ¹	15-20 miles
Turning radius.....	28"
Ground clearance	3.7"
Maximum step climb	2"
Length.....	47"
Width	27.5"
Drive tires	(2) 10" x 3.5"
Caster tires	(4) 8" x 2"
Total weight with batteries	273 lbs
Base only.....	114 lbs
Base with batteries.....	222 lbs
Uses 70 amp long-lasting batteries	

¹ Speed and range will vary with respect to terrain and weight of rider. The Manufacturer reserves the right to alter specifications without notice.

