



COBRA GT4 heavy-duty scooter



FEATURES & BENEFITS

- New, innovative design for a distinctive, stylish appearance.
- One-touch, infinite tiller angle adjustment allows the user to easily position the tiller into the ideal driving position.
- Ergonomic delta control handles makes it easier for users with limited strength and dexterity to operate the scooter.
- Large, air-filled, low profile tires and full suspension for smoothest outdoor ride quality.
- "Blade" style alloy wheel rims.
- Efficient, extra bright LED lighting. (Tail Lights & Turn Signals)
- Choice of Full or Sport Captain's Seats in 20" or 22" widths.
- Seats are height adjustable, swivel and feature a fold-down, reclining backrest. Seats can be adjusted forward or backward on auto-style slides.
- Padded armrests feature easy turn cylinders for angle adjustment.
- Positioning belt included.
- "High/Low" switch for instant speed control.
- Anti-tip wheels for greater safety.
- Easy, freewheel mode.

SPECIFICATIONS

Weight Capacity	450 lbs.	Batteries(2)	12V x 75AH
Top Speed	10 mph	Charger	5A Offboard
Max Range.....	35 miles	Brakes.....	Electromagnetic
Turning Radius.....	82"	Front Wheel.....	14" x 4" Flat-free
Climbing Angle.....	8°	Rear Wheels	14" x 4" Flat-free
Ground Clearance.....	7"	Base Weight	191 lbs.
Dimensions	63" L x 28.8" W	Battery Weight(2)	94 lbs.
Floor-to-Seat Height.....	27" - 30"	Seat Weight.....	55 lbs.
Seat Type.....	Captain's	Total Weight	340 lbs.
Seat Dimensions	20" W x 20" D / 22" W x 20" D	Warranty on Frame	Lifetime
Controller.....	Dynamic Rhino II 160A	Warranty on Electronic	
Motor.....	24V x 950W x 5100 rpm	Controller/Drive Train Components	14 Mths
		Warranty on Batteries	6 Mths