



TEKNIQUE® FWD

Front-wheel drive

The Hoveround Teknikue® FWD power chair is designed to easily maneuver both in the home and outdoors. Everyone will appreciate the comfortable design and powerful front-wheel drive maneuverability. The Teknikue FWD makes a great choice for improved mobility.

FEATURES & BENEFITS

- The better choice for individuals over 6 feet tall.
- Flip-up footplate and compact design for easy maneuverability.
- Higher ground clearance and powerful front-wheel drive performance.
- Adjustable width armrests for comfort and maximized support.
- Lightweight, portable charger for recharging power chair batteries anywhere.

SPECIFICATIONS

Max weight capacity	300 lbs
Max speed ¹	4.2 mph
Max range per charge ¹	12 miles
Turning radius	24"
Ground clearance	3.6"
Maximum step climb	2"
Length	41.9"
Width	25"
Drive tires	(2) 10" x 3.5"
Caster tires	(2) 8" x 2"
Total weight with batteries	176 lbs
Base only	88 lbs
Base with batteries	123 lbs

¹ Speed and range will vary with respect to terrain and weight of rider. The Manufacturer reserves the right to alter specifications without notice.