



The HOVERGlide HD stair lift by Hoveround has a 600 lb. capacity. High-backed, contoured 25" wide seat ensures a comfortable and secure ride. Sturdy armrests and reinforced footrest add stability for safer transfers. Slides out from wall as it pivots at top landing.

specifications

HOVERGlide HD

- Patented helical worm drive for safe, smooth, energy-efficient ride
- High capacity seat, arm and footrest
- Folding footrest
- Integrated constant-pressure controls
- Swivels at top landing for entry and exit
- Powerful enough to travel up to a 45-degree incline
- DC drive even runs during power interruptions
- Clean installation and maintenance; no grease or lubricants to apply or clean up
- 600 lb. weight capacity



The energy-efficient patented worm drive requires no greases or lubricants.



Safety sensors on carriage and platform stop the lift instantly at obstacles.

Specifications	SL600/SL600S
Maximum Weight Capacity	600 lbs. (272 kg)
Motor Output Amps at Rated Load	16 amps
Overspeed Governor	Yes
Track (rail) Type	Clear Anodized Aluminum
Travel	16' Standard; 40' Maximum (12.2 m)
Average Number of Return Trips per Charge	10 - 20
Symmetrical Lift	Yes
Control in Armrest	Yes
Need to Open Chassis to Install or Change Sides	No
Lift Mounts to Steps or Wall	Steps
Minimum Folded Width	17" (43 cm)
Minimum Footrest Height	3.6" (9 cm)
Clear Distance Between Armrests	25" (63.8 cm)
Floor to Seat Height Low / High	22"
Minimum Wall to Stair Side of Rail	12" (30.5 cm)
Seat Depth	16" (40.6 cm)
Backrest Height From Top of Seat	18" (45.7 cm)
Electrical Requirements	120VAC 15A
Operation Power	24V DC Battery
Speed	14 fpm
Incline Range	27 - 45°
Drive System	Helical Worm Gear
Safety Devices	Centrifugal Overspeed Brake Bi-directional Footrest obstruction Sensor Swivel Seat Cut-off Switch Final Limits Seat Belt Chassis Obstruction Sensors
Warranty:	3-Year on Components 10-Year on gear rack 1-Year on Batteries

Call today at 800-583-1204 or visit us online at www.hoveround.com