



TEKNIQUE® XHD

bariatric support 301-450 lbs

The highly maneuverable Teknik® XHD is the definitive bariatric power chair for anyone in need of a more robust mobility solution. The XHD is designed with a 450 lb weight capacity and runs on heavy-duty batteries for extended use. Combined with wider seating for comfort and tandem rear tires. The Teknik XHD stands out as the bariatric power chair of choice in its class.



FEATURES & BENEFITS

- Durable, heavy-duty construction for bariatric users from 301 - 450 lbs.
- Adjustable height and back positioning for extra comfort.
- Custom-contoured wide seat (22" or 24") for maximum riding support and comfort.
- Adjustable width armrests for comfort and proper upper body support.
- Oversized footplate.
- Lightweight, portable charger for recharging power chair batteries anywhere.
- Optional recline seat upgrade offers a higher seat back and adjustable recline providing additional comfort and support.*
- Optional head rest, oxygen tank, cane, crutches holder, elevated leg rests and many more accessories are available to customize your power chair.*

SPECIFICATIONS

Max weight capacity	450 lbs
Max speed ¹	4 mph
Max range per charge ¹	15 miles
Turning radius	28"
Ground clearance	3.7"
Maximum step climb	1.5"
Length	45"
Width	28.5"
Drive tires	(2) 10" x 3.5"
Caster tires	(4) 8" x 2"
Total weight with batteries	256 lbs
Base only	114 lbs
Base with batteries	192 lbs

* Optional accessories may not be medically necessary as determined by your doctor or covered by Medicare or third party reimbursement.

¹ Speed and range will vary with respect to terrain and weight of rider. The Manufacturer reserves the right to alter specifications without notice.