

HOVEROUND®

Personal Mobility Solutions



HOVEROUND LX-5®

One of the most maneuverable power wheelchairs available today, the Hoveround LX-5® will get you where you want to go. The unique, Round for a Reason® design enables easy navigation around tight corners and through narrow spaces without damaging walls or furniture. The LX-5® can safely function indoors as well as outdoors, offering power wheelchair users greater flexibility to maintain their freedom and independence.



Features & Benefits

- Round base and compact design provides exceptional maneuverability in tight spaces.
- Custom contoured innerspring seat provides all-day comfort.
- Adjustable width armrests for comfort and maximized support.
- Includes lightweight, portable charger for recharging power chair batteries anywhere.
- Available in Majestic Red or Vivacious Blue.
- Optional recline seat upgrade offers a higher seat back and adjustable recline providing additional comfort and support.

Specifications

Max weight capacity	300 lbs.
Max speed ¹	4 mph
Max range per charge ¹	12 miles
Turning radius	22.7"
Ground clearance	2.5"
Maximum step climb	2"
Length	38"
Width	24"
Drive tires	9" x 3.5"
Caster tires	8" x 2"
Total weight with batteries	173 lbs.
Base only	88 lbs.
Base with batteries	123 lbs.

¹ Speed and range will vary with respect to terrain and weight of rider. The Manufacturer reserves the right to alter specifications without notice.

Don't let a life with limited mobility limit the life you live!