



HOVEROUND MPV5®

Most maneuverable

The MPV5® is not just one of the most maneuverable power chairs available today, it is also Hoveround's most popular model. The unique, Round for a Reason® design enables easy navigation around tight corners and through narrow spaces without damaging walls or furniture. The MPV5 can safely function indoors as well as outdoors, offering power wheelchair users greater flexibility in their freedom and mobility.

FEATURES & BENEFITS

- Flip-up footplate, round frame base and compact design for easy maneuverability.
- Custom-contoured inner-spring seat provides for all day comfort.
- Adjustable width armrests for comfort and maximized support.
- Lightweight, portable charger for recharging power chair batteries anywhere.
- Shown with optional red color panels which are available for purchase.†

SPECIFICATIONS

Max weight capacity	300 lbs
Max speed¹	4 mph
Max range per charge¹	12 miles
Turning radius	22.7"
Ground clearance	2.5"
Maximum step climb	2"
Length	38"
Width	24"
Drive tires	(2) 9" x 3.5"
Caster tires	(2) 8" x 2"
Total weight with batteries	173 lbs
Base only	88 lbs
Base with batteries	123 lbs

† Optional color panels are not available for Medicare or third party reimbursement. May be purchased separately.

¹ Speed and range will vary with respect to terrain and weight of rider. The Manufacturer reserves the right to alter specifications without notice.