




PHOENIX HD 4WHL heavy-duty travel scooter



FEATURES & BENEFITS

- Quick and easy disassembly.
- Includes interchangeable panels in Red & Blue.
- Armrests are padded and adjustable.
- Quick connect batteries.
- Ergonomic throttle control.
- Height-adjustable swivel seat with two-tone upholstery.

- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Headlight.
- Swivel seat with fold-down backrest.
- Easy freewheel operation.

SPECIFICATIONS

HCPCS.....	K0801
Weight Capacity	350 lbs.
Top Speed.....	4 mph
Max Range.....	15 miles
Turning Radius.....	54"
Climbing Angle.....	6°
Ground Clearance.....	3.75"
Dimensions	41.5" L x 22.5" W
Floor-to-Seat Height.....	20" - 22"
Seat Type.....	Folding
Seat Dimensions	17.5" W x 17" D
Motor.....	24V x 350W x 3600 rpm
Batteries(2)	12V x 20AH

Charger	1.5A Offboard
Brakes.....	Electromagnetic
Front Wheel.....	8" x 2" Flat-free
Rear Wheels	8" x 2" Flat-free
Base Weight	77 lbs.
Battery Weight(2)	31 lbs.
Seat Weight.....	21 lbs.
Total Weight	129 lbs.
Heaviest Piece Weight.....	46 lbs. (Front section)
Warranty on Frame	Lifetime
Warranty on Electronic	
Controller/Drive Train Components	14 Mths
Warranty on Batteries	6 Mths