



## SCOUT 3WHL

### compact travel scooter



### FEATURES & BENEFITS

- Quick and easy disassembly.
- Interchangeable red and blue panels.
- Armrests are padded and adjustable.
- Quick connect batteries.
- Ergonomic throttle control.
- Height-adjustable swivel seat with fold-down backrest.

- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Delta tiller.
- Easy freewheel operation.

### SPECIFICATIONS

HCPCS.....	K0800
Weight Capacity .....	300 lbs.
Top Speed .....	4.25 mph
Max Range.....	9 miles
Turning Radius.....	45.5"
Climbing Angle.....	6°
Ground Clearance.....	2.5"
Dimensions .....	42.5" L x 20.5" W
Floor-to-Seat Height .....	17" - 19"
Seat Type.....	Folding
Seat Dimensions .....	16.5" W x 16" D
Controller.....	CTE 45A
Motor.....	24V x 270W x 4700 rpm
Batteries(2) .....	12V x 12AH

Charger .....	2A Offboard
Brakes.....	Electromagnetic
Front Wheel.....	8" x 2" Flat-free
Rear Wheels .....	8" x 2" Flat-free
Base Weight .....	49 lbs.
Battery Weight(2) .....	20 lbs.
Seat Weight.....	16 lbs.
Total Weight .....	85 lbs.
Heaviest Piece Weight.....	29 lbs.
Warranty on Frame .....	Lifetime
Warranty on Electronic	
Controller/Drive Train Components .....	14 Mths
Warranty on Batteries .....	6 Mths