



TEKNIQUE RWD

Rear-wheel drive

The Teknique® RWD custom power wheelchair is ideal for use over many types of terrain. Featuring the added support of swing-away footrests and rear wheel drive mobility, the Teknique RWD is the perfect power wheelchair for anyone planning to be out and about most of the day.

FEATURES & BENEFITS

- Fixed-angle, swing-away footrests with heel loops for secure leg positioning.
- Sturdy cup holder and side tote bag for increased independence and convenience.
- Adjustable armrests for comfort and proper upper body support.
- Custom-contoured seat for maximum comfort.
- Dual motors with powerful rear-wheel drive for extra riding stability.
- Lightweight, portable charger for recharging power chair batteries anywhere.

SPECIFICATIONS

Max weight capacity.....	300 lbs
Max speed ¹	5.5 mph
Max range per charge ¹	15-20 miles
Turning radius.....	36"
Ground clearance.....	2.5"
Maximum step climb.....	2"
Length.....	42"
Width.....	25"
Drive tires.....	(2) 10" x 3.5"
Caster tires.....	(2) 8" x 2"
Total weight with batteries.....	229 lbs
Base only.....	100 lbs
Base with batteries.....	178 lbs

¹ Speed and range will vary with respect to terrain and weight of rider.
The Manufacturer reserves the right to alter specifications without notice.