



TRANSPORTER GL

Sporty & Stylish

Keep your active lifestyle with our sporty and stylish Transporter® GL scooter. Designed primarily for those in search of a little added assistance with mobility, this 3-wheel scooter is perfect for getting you around town safely. Affordable, and equipped with a 360-degree swivel, sliding seat, basket, and adjustable tiller, the Transporter GL is one of the most user-friendly electric scooters you'll find on the market today.

FEATURES & BENEFITS

- Rugged, three-wheel, compact design with front light makes it sporty, stylish and safe.
- 360-degree swivel, sliding seat for maximum convenience.
- Adjustable back angle with ergonomic captain's seat and headrest for utmost comfort.
- Steering tiller adjusts for a more individually customized ride.
- Lightweight, portable charger lets you charge your batteries anywhere.
- Attachable basket for transporting items while on-the-go.

SPECIFICATIONS

Max weight capacity.....	300 lbs
Max speed ¹	4 mph
Max range per charge ¹	25 miles
Turning radius.....	39"
Ground clearance.....	3"
Length.....	45"
Width.....	24"
Front tires.....	(1) 9" x 3"
Back tires.....	(2) 11" x 3.5"
Total weight with batteries.....	178 lbs

¹ Speed and range will vary with respect to terrain and weight of rider.
The Manufacturer reserves the right to alter specifications without notice.