



TRANSPORTER

Sporty & Compact

The sporty Transporter® is one of the most compact mobility scooters for getting about and enjoying your independence. This affordable and easy-to-manuever 4-wheel scooter has several convenient features that add comfort, stability and safety to your ride. Now, do more than just imagine a sunny day at the park or dream of visiting your neighbors down the street. With the Transporter mobility scooter, you have the power to take yourself there!

FEATURES & BENEFITS

- Adjustable tiller with ergonomic hand grips for ease of operation.
- Easy-to-use key start with speed adjustments.
- Padded, 17" swivel seat and flip-up armrests for added comfort and safer access.
- Durable, flat-free, non-marking tires and anti-tip wheels for increased stability.
- Convenient plastic basket with removable liner for carrying personal items.
- Disassembles without tools and back rest folds down for convenient transport.
- Includes set of three interchangeable body panels in red, blue and silver so you can customize your scooter.
- 15 miles maximum range with optional 21 AH battery pack.

SPECIFICATIONS

Max weight capacity.....	300 lbs
Max speed ¹	4 mph
Max range per charge ¹	9 miles w/ 12AH
Turning radius.....	43"
Ground clearance	2"
Maximum step climb	2"
Length.....	42"
Width	19"
Front tires.....	(2) 7" x 2"
Back tires.....	(2) 8" x 2"
Total weight with batteries	103 lbs

¹ Speed and range will vary with respect to terrain and weight of rider.
The Manufacturer reserves the right to alter specifications without notice.