



**VENTURA 3WHL**  
midsize scooter



**FEATURES & BENEFITS**

- Includes interchangeable color panels in Red & Blue.
- Armrests are padded and adjustable.
- Floorboard with 23" of foot space.
- Ergonomic throttle control.
- Height-adjustable swivel seat with fold-down backrest.
- Flat-free, non-marking tires.
- Seat position can be adjusted forward or backward.
- Rear suspension.
- Large, plastic carry basket.
- Easy-to-adjust Delta tiller.
- Anti-tip wheels.
- Headlight.
- Easy freewheel operation.

**SPECIFICATIONS**

Weight Capacity.....	350 lbs.	Charger.....	4A Offboard
Top Speed.....	5 mph	Brakes.....	Electromagnetic
Max Range.....	15 miles	Front Wheel.....	10" x 3" Flat-free
Turning Radius.....	53"	Rear Wheels.....	10" x 3" Flat-free
Climbing Angle.....	7.5°	Base Weight.....	88 lbs.
Ground Clearance.....	2.5"	Battery Weight(2).....	28 lbs.
Dimensions.....	48" L x 25" W	Seat Weight.....	26 lbs.
Floor-to-Seat Height.....	21.5" - 23.5"	Total Weight.....	142 lbs.
Seat Type.....	Folding	Warranty on Frame.....	Lifetime
Seat Dimensions.....	18" W x 17.5" D	Warranty on Electronic	
Motor.....	24V x 350W x 5100 rpm	Controller/Drive Train Components.....	14 Mths
Batteries(2).....	12V x 20AH	Warranty on Batteries.....	6 Mths