



VENTURA 4WHL
midsize scooter



FEATURES & BENEFITS

- Includes interchangeable panels in Red & Blue.
- Armrests are padded and adjustable.
- Ergonomic throttle control.
- Height-adjustable swivel seat with fold-down backrest.
- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Seat position can be adjusted forward or backward.

- Front and rear suspension.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Headlight.
- Delta tiller.
- Easy freewheel operation.

SPECIFICATIONS

Weight Capacity.....350 lbs.
 Top Speed..... 4 mph
 Max Range.....15 miles
 Turning Radius.....54"
 Climbing Angle..... 7.5°
 Ground Clearance.....2.5"
 Dimensions.....48" L x 25" W
 Floor-to-Seat Height.....22" - 24"
 Seat Type.....Folding
 Seat Dimensions.....18" W x 17.5" D
 Motor.....24V x 350W x 5100 rpm
 Batteries(2).....12V x 20AH

Charger.....4A Offboard
 Brakes..... Electromagnetic
 Front Wheel..... 10" x 3" Flat-free
 Rear Wheels..... 10" x 3" Flat-free
 Base Weight.....99 lbs.
 Battery Weight(2).....28 lbs.
 Seat Weight.....26 lbs.
 Total Weight.....153 lbs.
 Warranty on Frame..... Lifetime
 Warranty on Electronic
 Controller/Drive Train Components..... 14 Mths
 Warranty on Batteries..... 6 Mths