



VENTURA DLX 3WHL midsize scooter



FEATURES & BENEFITS

- Includes interchangeable color panels in Red & Blue.
- 36AH batteries for extended range.
- Armrests are padded and adjustable.
- Floorboard with 23" of foot space.
- Ergonomic throttle control.
- Height-adjustable swivel seat features fold-down backrest.
- Seat position can be adjusted forward or backward.
- Rear suspension.
- Programmable controller.

- Seat back reclines and features adjustable headrest.
- Available with 18" or 20" width Captain's Seat.
- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Easy-to-adjust Delta tiller.
- Anti-tip wheels.
- Headlight.
- Rearview mirrors.
- Turn signals/hazard lights.
- Easy freewheel operation.

SPECIFICATIONS

HCPCS.....	N/A
Weight Capacity	350 lbs.
Top Speed	5 mph
Max Range.....	22 miles
Turning Radius.....	.53"
Climbing Angle.....	7.5°
Ground Clearance.....	2.5"
Dimensions	48" L x 25" W
Floor-to-Seat Height.....	20.5" - 23.5"
Seat Type.....	Captain's
Seat Dimensions	18" W x 18" D / 20" W x 18" D
Motor.....	24V x 350W x 5100 rpm
Batteries(2)	12V x 36AH

Charger	4A Offboard
Brakes.....	Electromagnetic
Front Wheel.....	10" x 3" Flat-free
Rear Wheels	10" x 3" Flat-free
Base Weight	88 lbs.
Battery Weight(2)	46 lbs.
Seat Weight.....	40 lbs.
Total Weight	174 lbs.
Warranty on Frame	Lifetime
Warranty on Electronic Controller/Drive Train Components	14 Mths
Warranty on Batteries	6 Mths