



TEKNIQUE® XHD

Bariatric support

The highly maneuverable Teknique® XHD is the definitive bariatric power chair for anyone in need of a more robust mobility solution. The XHD is designed with a 450-pound weight capacity and runs on heavy-duty batteries for extended use. Combined with wider seating for comfort and tandem rear tires. The Teknique XHD stands out as the bariatric power chair of choice in its class.

FEATURES & BENEFITS

- Durable, heavy-duty construction for bariatric users from 301 - 450 lbs.
- Adjustable height and back positioning for extra comfort.
- Custom-contoured wide seat (22" or 24") for maximum riding support and comfort.
- Adjustable width armrests for comfort and proper upper body support.
- Oversized footplate.
- Lightweight, portable charger for recharging power chair batteries anywhere.

SPECIFICATIONS

Max weight capacity	450 lbs
Max speed ¹	4 mph
Max range per charge ¹	15 miles
Turning radius	28"
Ground clearance	3.7"
Maximum step climb	1.5"
Length	45"
Width	28.5"
Drive tires	(2) 10" x 3.5"
Caster tires	(4) 8" x 2"
Total weight with batteries	256 lbs
Base only	114 lbs
Base with batteries	192 lbs

¹ Speed and range will vary with respect to terrain and weight of rider. The Manufacturer reserves the right to alter specifications without notice.