

AEROLUX 2Comfort

ULTRALIGHT, ULTRA COMFORTABLE



**100%
CARBON FIBRE**



**HOLDS UP TO
300 LBS**




**WEIGHS ONLY
33 LBS**

EXPERIENCE NEXT-LEVEL COMFORT

- Ergonomic padded backrest delivers all-day comfort.
- Larger seating area, with memory foam waterfall cushion (18"W x 17"D).
- Extra wide padded armrests offer enhanced support and cushioning.
- Market leading 18.5" wide between armrests, and extendable to 19.5" with spacers (included).
- 3-position adjustable footrest reduces fatigue.
- Larger wide-stance wheels deliver the best ride quality and stability in any carbon chair.
- Ultra quiet 200W motors and power delivery.
- Same unique cross-folding frame as the original AEROLUX.
- Optional Bluetooth remote control with extended range.



KEY SPECIFICATIONS



ACCESSORIZE YOUR AEROLUX 2C

Enhance your ride with coloured seating sets and various add-ons!



PURPLE



RED



BLUE

