

merits

# ATLANTIS P710A

Super Heavy Duty RWD Powerbase

K0826 / K0827



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- Heavy-duty power wheelchair with a weight capacity of 600 lbs.
- Four seat posts provide a stable and durable structure for the frame
- Extra-wide seat and extra-high back provide great comfort for bariatric users
- Super high power in-line motors help riders negotiate 10° uphill ramps with ease

## Accessories

- Oxygen tank holder
- Cane & crutch holder
- Walker holder
- Foldable rear basket
- Cup holder
- Side bag
- USB charger
- Cover

## Options

- Captain power seat option
- Elevating legrest available (K0195)
- Swing-away joystick
- Flip-up, height-adjustable, detachable armrests



## Specifications

Product#	P710A
Weight Capacity	600 lbs.
Drivetrain	RWD
Seat Type	Captain
Seat Width	24.5"
Seat Depth	20"
Backrest Height	21.25"
Seat-to-Floor Height	22" - 23"
Overall Length	33"
Overall Width	32"
Overall Height	41.5"
Turning Radius	31.5"
Ground Clearance	2.75"
Total Weight	304 lbs.
Wt of Heaviest Pc.	167 lbs.
Max Speed	5 mph
Range Per Charge	up to 32 mi.
Motor	DC 24V / 300W
Controller	Dynamic LiNX
Battery Size (2 required)	75Ah, M24 AGM
Charger	6A off-board
Gradient	10°
Front Wheels	9" pneumatic
Drive Wheels	12" pneumatic
Rear Wheels	3" anti-tippers
HCPCS Codes	K0826 / K0827
Color Options	<div style="display: inline-block; width: 15px; height: 15px; background-color: red; margin-right: 5px;"></div> Candy Apple Red <div style="display: inline-block; width: 15px; height: 15px; background-color: blue; margin-left: 10px; margin-right: 5px;"></div> Metallic Blue

Features MK batteries



merits



**merits** Enhancing life through mobility.

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0425