

Cobra GT4 Heavy Duty Power Mobility Scooter

HCPCS:



From its distinctive look to class-leading performance, we have stopped at nothing to create the new standard in heavy duty scooters. The Cobra GT4 boasts a dynamic, modern design and incorporates cutting edge technological advancements to provide the most attractive, highest performance scooter available. In addition, the Cobra provides maximum comfort with a smooth, stable ride and a deluxe Captain's Seat with a width of 22". Experience a new age in scooter design with the Cobra GT4.

- Ergonomic delta control handles makes it easier for users with limited strength and dexterity to operate the scooter
- Innovative design for a distinctive, stylish appearance
- One-touch, infinite tiller angle adjustment allows the user to easily position the tiller into the ideal driving position
- Large, air-filled, low profile tires and full suspension for smoothest outdoor ride quality
- "Blade" style alloy wheel rims
- Efficient, extra bright LED lighting
- Captain's seat is height adjustable, swivel and features a fold-down, reclining backrest; seat can be adjusted forward or backward on auto-style slides
- Padded armrests feature easy turn cylinders for angle adjustment
- Positioning belt included
- "High/Low" switch for instant speed control
- Anti-tip wheels for greater safety
- Easy, freewheel mode

PRODUCTS	Description
COBRAGT422CS	COBRA GT4 SCOOTER WITH 22" CAPTAIN SEAT
SPECIFICATIONS	Description
Warranty	24 Month Limited Warranty on Electronics
Product Weight	340 lb
Product Dimensions	63" (L) x 29" (W) x 53" (H)
Seat Dimensions	20" (D) x 22" (W)
Weight Capacity	450 lb
Top Speed	10 mph
Batteries (pair)	12V x 75AH
Battery Charger	5A Offboard
Climbing Angle	8°
Color	Silver
Controller	Dynamic Rhino 2 120A
Ground Clearance	5.5"
HCPCS	N/A

SPECIFICATIONS	Description
Maximum Range	22 miles
Motor	24V x 950W x 5200 RPM
Seat Type	Captain - Sport Style
Turning Radius	82"
Warranty on Batteries	12 Month Warranty on Batteries
Warranty on Frame	Limited Lifetime Warranty on Frame
Wheel Size	14"