

Model	Go Chair® MED
Weight Capacity	300 lbs.
Maximum Speed <sup>2</sup>	Up to 3.7 mph
Ground Clearance <sup>4</sup>	1.75" at motor
Turning Radius <sup>4</sup>	25.5"
Overall Length <sup>4</sup>	33" with footplate
Overall Width <sup>4</sup>	25.5" with armrests
Seat-to-Ground Height	22"-23"
Front Wheels	6" solid casters
Drive Wheels	9" solid
Rear Wheels	2" solid anti-tips
Suspension	Limited
Range Per Charge <sup>2,3</sup>	Up to 9.41 miles
Heaviest Piece When Disassembled <sup>4</sup>	36 lbs. (rear section)
Standard Seating <sup>6</sup>	Type: Compact, foldable Material: Black vinyl Dimensions: Width: 18", Depth: 17", Weight: 29 lbs. or Width: 20", Depth: 18", Weight: 32 lbs.
Drivetrain	Two-motor, rear-wheel drive
Braking System	Regenerative and electromechanical
Standard Electronics	40 amp, PG GC3
Battery Requirements <sup>5</sup>	20 Ah battery pack Weight: 31.4 lbs.
Battery Charger	Off-board, 2 amp
Warranty	Frame: 5-year limited Drivetrain: 13-month limited Electronics: 13-month limited Batteries: 13-month limited

### HCPCS Code: Group 2: Standard, Portable (K0821)

- 1 Pride FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 3 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 5 Battery weight may vary based on battery manufacturer.
- 6 Seat weight includes the joystick.

**PRIDE  
MOBILITY®**  
800-800-8586 • [pridemobility.com](http://pridemobility.com)



HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

082525

**GOCHAIR®  
MED**

**TRAVEL  
MADE SIMPLE**



FDA Class II  
Medical Device¹

800-800-8586 • [pridemobility.com](http://pridemobility.com)



# POWERFUL, PORTABLE INDEPENDENCE.




FDA Class II Medical Device<sup>1</sup>



Disassembles quickly for  
transport or storage



Optional swing-away joystick  
swings out at an angle for  
convenient access



Available with convenient  
cup holder

The Go Chair® MED is #1 for portability! With convenient feather-touch disassembly, the Go Chair MED allows you to experience lightweight travel and independence on the go. Compact, portable and perfect for travel, the Go Chair® MED is a feature-rich power chair that disassembles easily in just a few steps!

## FEATURES

- 18" x 17" or 20" x 18", 5-point swivel seat with latch release and 60" lap belt
- Compact turning radius
- Large perforated footplate
- One-hand feather-touch disassembly
- The heaviest piece is only 36 lbs.

## OPTIONS

- Swing-away joystick for convenient access
- Medium back memory seat options available in 18" x 18" or 20" x 18"

## AVAILABLE COLORS

Black

## OPTIONAL ACCESSORIES\*

- Cane crutch holder (single/double rear mount)
- Cup holder
- Oxygen tank holder
- Lap belts in 50" or 70"
- Saddle bag
- Safety flag
- USB charger
- Walker holder
- Weather cover
- Military Patches

*\* Only one rear accessory can be used at a time.*