

**G**

## GRAND TOURER

MOBILITY HAS NEVER LOOKED THIS GOOD



**Users Can Easily Lift  
Into Cars With Built-in  
Loading Wheels**

**3 Distinct Driving  
Modes That Adapt  
to Any Setting**

**5 Adjustable  
Armrest Height  
Levels**

Designed to turn heads, the Grand Tourer (GT) combines standout style with uncompromising performance. Thoughtfully engineered to address the most common user pain points, the GT features exceptional comfort, smooth ride quality, effortless loading, and a modern, sleek style. Because yes, you really can have it all.

# STYLE YOU'LL LOVE, PERFORMANCE YOU'LL TRUST.

- **Power Meets Precision** — With 250W Sine Wave Motors, enjoy a ride that's powerful, quiet, and incredibly responsive.
- **Advanced Sine Wave Controller** — A sleek LED display with intuitive vibration feedback and effortless mode switching delivers smoother handling and a more responsive ride.
- **3 Distinct Driving Modes** — With Soft, Normal, and Sport modes, the GT adapts effortlessly to any setting, whether it's gentle control indoors, balanced daily performance, or maximum agility outdoors.
- **Effortless Loading** — Loading wheels that lock at 90° or 135° allow you to simply tilt and push the chair into any vehicle. This significantly reduces the chair weight, making loading possible for every user.
- **Quick & Secure Folding** — The quick-fold handle makes storage simple, and the safety lock keeps everything secure until you're ready to unfold.
- **Adjustable Height** — With 5 armrest settings, a backrest that extends up to 4" higher, and a footrest with 2.4" of added range, every ride feels precisely tailored.
- **Seamless Access** — Fold-up armrests create a fully unobstructed side entry, ensuring easier transfers and improved accessibility.
- **USB-C Charging Port** — Stay powered on the go with the built-in USB-C port in the joystick, perfect for charging your phone or other devices while you ride.



## KEY SPECIFICATIONS

<b>CHAIR NET WEIGHT</b>	49 LBS	<b>INCLINE CAPABILITY</b>	12°	<b>FOLDED DIMENSIONS</b>	29.5"L X 23.6"W X 13.8"H
<b>CHAIR GROSS WEIGHT</b>	53 LBS	<b>TURNING RADIUS</b>	25"	<b>UNFOLDED DIMENSIONS</b>	39"L X 23.6"W X 40"H
<b>MAX WEIGHT CAPACITY</b>	300 LBS	<b>MOTOR POWER</b>	2X 250 W BRUSHLESS	<b>SEAT CUSHION SIZE</b>	18"W X 18"D X 3"H
<b>MAX DRIVING RANGE</b>	15.5 MI / 25 KM	<b>BATTERY TYPE</b>	Li-Ion 24 V / 12 Ah	<b>FRONT WHEEL DIAMETER</b>	8"
<b>MAX SPEED</b>	3.7 MPH / 6 KPH	<b>SEAT TO FLOOR HEIGHT</b>	20"	<b>REAR WHEEL DIAMETER</b>	10"

## GT ACCESSORIES



STORAGE BAG



BATTERY &  
JOYSTICK BAG



CUP HOLDER

## CHOOSE FROM 2 SLEEK FRAME COLOURS



MATTE BLACK



MATTE WHITE