



ULTRA CST P325

MWD Powerbase | Captain Seat with Tilt

K0836



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features


- Center of gravity tilt 0° - 50° posterior
- Premium, semi-reclining ergonomic captain seat
- Tension-adjustable front and rear suspension
- Adjustable seat width and depth from 18" - 20"; depth adjustable in 1" increments
- 24" base width and seat-to-floor height as low as 21.5"
- Rated up to 19 miles per charge

Options

- Power articulating foot platform
- Elevating leg rests or swing-away footrests
- Swing-away joystick
- USB charger



Specifications

Product#	P325
Weight Capacity	300 lbs.
Drivetrain	MWD
Maximum Speed	5 mph
Per-Charge Range	11 mi. / 19 mi.
Battery Size (2 required)	55Ah, 22NF AGM 51Ah, 22NF GEL 33Ah, U1 AGM
Ground Clearance	2.5"
Turning Radius	28"
Overall Length	36"
Overall Width	24"
Overall Height	40"
Suspension	Front and Rear
Seat-to-Floor Height	21.5" - 22.5" (no cushion)
Seat Width	18"-20", adjusted in 1" increments
Seat Depth	18"-20", adjusted in 1" increments
Backrest Height	20"
Total Weight	310 lbs.
Motor	DC 24V, 240V
Controller	P&G VR2-A2 60A
Charger	6A off board
Gradient	6°
Front Wheels	6" x 2" solid
Drive Wheels	14" x 3" foam-filled
Casters	6" x 2" solid
HCPCS Codes	K0836
Color Options	 Red

Features MK batteries



 **Enhancing life through mobility.**

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0425