



ULTRA P325

Captain seat | Pan seat | Single Power

K0822 / K0823 / K0835



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Center of gravity tilt 0° - 50° posterior
- Shear reduction recline 90° - 110°
- Tension-adjustable front and rear suspension
- Adjustable seat width and depth from 16" - 22"; depth adjustable in 1" increments
- 24" base width and seat-to-floor height as low as 16.5"
- Rated up to 19 miles per charge

Options

- 12" elevating seat
- Power articulating foot platform
- Elevating leg rests or swing-away footrests
- Swing-away joystick
- Oxygen tank holder
- USB charger

Models

Vision Ultra - G2, Standard, Captain P325-2S C

Vision Ultra - G2, Standard, Pan P325- 2S S

Vision Ultra - G2, Standard, Single Power P325-2SSP S



Specifications

Product#	P325
Weight Capacity	300 lbs.
Drivetrain	MWD
Maximum Speed	5 mph
Per-Charge Range	12 mi. / 19 mi.
Battery Size (2 required)	55Ah, 22NF AGM 51Ah, 22NF GEL 33Ah, U1 AGM
Ground Clearance	2.5"
Turning Radius	26"
Overall Length	36"
Overall Width	24"
Overall Height	40"
Suspension	Front and Rear
Seat-to-Floor Height	Power tilt or pan: 16.5"-19.5", Captain: 20"-23"
Seat Width	16"-22", adjusted in 1" increments
Seat Depth	16"-22", adjusted in 1" increments
Backrest Height	18"-28"
Total Weight	315 - 384 lbs.
Motor	2-Pole/240W
Controller	P&G VR2-A2 60A
Charger	6 A off board
Gradient	6°
Front Wheels	6" x 2" foam-filled
Drive Wheels	14"x 3" foam-filled
Casters	6" x 2" foam-filled
HCPCS Code	K0822 / K0823 / K0835
Color Options	<input type="checkbox"/> White <input checked="" type="checkbox"/> Red <input checked="" type="checkbox"/> Blue <input checked="" type="checkbox"/> Charcoal Gray

Features MK batteries



merits Enhancing life through mobility.

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0425