

## Ventura 3-Wheel Scooter

HCPCS:



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

The Ventura 3 Wheel Scooter by Drive Medical comes with a Captain's Seat for superior comfort.

- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue (Figure A)
- 33 AH batteries for extended range
- Floorboard with 23" of foot space
- Ergonomic throttle control
- Height-adjustable swivel seat features fold-down backrest and adjustable padded armrests
- Seat position can be adjusted forward or backward
- Rear suspension and Anti-tip wheels
- Programmable controller
- Seat back reclines and features adjustable headrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust Delta tiller
- Headlight, Turn signals/hazard lights and Rearview mirrors
- Easy free wheel operation

PRODUCTS	Description
VENTURA320CS	VENTURA 3-WHL DELUXE 20"CAPT.SEAT,RED/BL
SPECIFICATIONS	Description
Warranty	24 Month Limited Warranty on Electronics
Product Weight	187 lb
Product Dimensions	48" (L) x 25" (W) x 45" (H)
Seat Dimensions	18" (D) x 20" (W)
Weight Capacity	400 lb
Brakes	Electromagnetic
Top Speed	6 mph
Batteries (pair)	12V x 33AH
Battery Charger	4A Offboard
Battery Weight (pair)	46 lbs.
Charger	4A Offboard
Climbing Angle	7.5°
Color	Red & Blue
Controller	PG 70A

SPECIFICATIONS	Description
Dimensions	48" (L) x 25" (W)
Ground Clearance	2.5"
HCPCS	N/A
Heaviest Piece Weight	60 lbs.
Maximum Range	16 miles
Motor	24V x 350W x 5100 RPM
Seat Type	Captain - Reclining High
Turning Radius	54"
Warranty on Batteries	12 Month Warranty on Batteries
Warranty on Frame	Limited Lifetime Warranty on Frame
Wheel Size	Front: 10" x 3"   Rear: 10" x 3"
ACCESSORIES	Description
SF8030	Crutch Holder Set, SType, 1/ea, EA