

|   |  |
|---|--|
| Model   | Victory® 10 <sup>1</sup>   |
| Product Number                                | 3W: SC610 4W: SC710  |
| Weight Capacity                               | 400 lbs.   |
| Maximum Speed <sup>2</sup>                    | Up to 5.3 mph  |
| Ground Clearance <sup>4</sup>                 | 3W: 2.2" at motor, 1.26" at anti-tip<br>3.4" at side deck<br>4W: 2" at motor, 1.5" at anti-tip<br>3.1" at side deck  |
| Turning Radius <sup>4</sup>                   | 3W: 42" 4W: 45.5"  |
| Overall Length <sup>4</sup>                   | 3W: 45.25" 4W: 47"   |
| Overall Width <sup>4</sup>                    | 3W: 22.25" 4W: 22.25"  |
| Seat-to-Ground Height                         | 23.5"-25"  |
| Seat-to-Deck Height                           | 17.5"-19"  |
| Front Tires                                   | 3.5" x 10" solid   |
| Rear Tires                                    | 3.5" x 10" solid   |
| Suspension                                    | N/A  |
| Range Per Charge <sup>2,3</sup>               | Up to 15.5 miles   |
| Total Weight <sup>6</sup>                     | Without batteries:<br>3W 133 lbs.; 4W: 137 lbs.<br>With batteries:<br>3W-181 lbs. (U1);<br>3W-197 lbs. (40AH);<br>4W-185 lbs. (U1);<br>4W-201 lbs. (40AH); |
| Heaviest Piece When Disassembled <sup>4</sup> | 3W: 54 lbs. (front section)<br>4W: 63 lbs. (front section)   |
| Standard Seating                              | Type: Compact, foldable<br>Weight: 28 lbs.<br>Material: Black Vinyl<br>Dimensions: 18" width, 17" depth  |
| Drive System                                  | Rear-wheel drive<br>sealed transaxle with 24-volt<br>DC motor  |
| Dual Braking System                           | Regenerative and<br>electromechanical  |
| Battery Requirements <sup>5,7</sup>           | Size: (2) U1 standard,<br>Weight: 24.5 lbs. each<br>Size: (2) 40 AH (optional),<br>Weight: 32.5 lbs. each  |
| Battery Charger                               | 3.5-amp, off-board (U1),<br>5-amp, off-board (40 AH)   |
| Warranty                                      | Frame: lifetime limited<br>Drivetrain: 2-year limited<br>Electronics: 2-year limited<br>Battery: 6-months  |

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

## PRIDE MOBILITY®

800-800-8586 • [pridemobility.com](http://pridemobility.com)



### HCPCS Code: K0801 (Group 1, Heavy Duty)

- <sup>1</sup> Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- <sup>2</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- <sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- <sup>4</sup> AGM or gel cell type required.
- <sup>5</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- <sup>6</sup> Battery weight may vary based on manufacturer.
- <sup>7</sup> Includes standard seat weight.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

081825



# COMFORT AND PERFORMANCE, NOW WITH TIGHTER TURNING!

FEATURING EXCLUSIVE PATENTED **EZ TURN**  
TECHNOLOGY



FDA Class II  
Medical Device<sup>1</sup>

800-800-8586 • [pridemobility.com](http://pridemobility.com)



# NOW WITH TIGHTER TURNING!

FEATURING EXCLUSIVE PATENTED



18" x 17"  
Pinchless hinge,  
black vinyl seat  
with 8-point  
swivel and  
sliders for  
increased  
comfort



Ambient  
voltmeter

Standard  
front basket

Easy-to-grip  
tiller adjustment  
knob

Long-lasting,  
maximum-  
intensity LED  
headlight  
offers optimal  
pathway  
illumination

Exclusive stylish, lightweight,  
non-scuffing, black, low-profile tires

The **Victory® 10** series offers added style, comfort & performance for your active lifestyle. Standard features include LED pathway light, LED battery meter, feather touch disassembly, wraparound delta tiller, solid tires and a front basket. A 400 lb. weight capacity, a per charge range of up to 15.5 miles and a maximum speed up to 5.3 mph compliment the versatility of the Victory 10.

## FEATURES

- Pride's exclusive EZ Turn technology (4W unit only) provides **3-wheel turning ability** while maintaining **4-wheel stability**
- Exclusive one-hand feather touch disassembly
- Charger port located on the tiller lets you conveniently charge your scooter
- Delta tiller with ergonomic wraparound handles lets you operate the scooter with one hand and rest your wrist
- Easy-to-read ambient voltmeter
- Long-lasting, maximum-intensity LED headlight offers optimal pathway illumination

## OPTIONAL SEATING

Limited recline, black vinyl, high back with height adjustable headrest, flip-up height and width adjustable armrests

- 18" x 18" (sliders) • 20" x 20" • 22" x 18"

Deluxe contour high back black vinyl with height adjustable headrest, flip-up height and width adjustable armrests

- 18" x 18" • 22" x 20"

Pinchless black vinyl medium back hinge seat with 8-point swivel and sliders for increased comfort

- 20" x 18"

## AVAILABLE COLORS

Candy Apple Red

Viper Blue

## OPTIONAL ACCESSORIES\*

- |                                       |  |
|---------------------------------------|--|
| • Cane clip attachment                | • Safety flag  |
| • Cane/crutch holder (single, double) | • Walker holder  |
| • Cell phone holder                   | • Weather cover  |
| • Cup holder                          | • XLR USB charger  |
| • Oxygen tank holder                  | • Lap belts: 50", 60" and 70"  |
| • Rear basket                         | <b>Cup holder, rear basket, and saddlebag can be ordered as a bundle</b> |
| • Rearview mirror(s)                  |  |
| • Saddle bag                          |  |

\* Only one rear accessory can be used at a time.



Tiller console



Feather-touch disassembly



Exclusive EZ Turn Technology  
on 4-Wheel Model

4-wheel model shown;  
3-wheel model available