



TEKNIQUE® FWD

comfort and support

The Hoveround Teknique® FWD power chair is designed to easily maneuver both in the home and outdoors. Everyone will appreciate the comfortable design and powerful front-wheel drive maneuverability. The Teknique FWD makes a great choice for improved mobility.



FEATURES & BENEFITS

- The better choice for individuals over 6' tall.
- Flip-up footplate and compact design for easy maneuverability.
- Higher ground clearance and powerful front-wheel drive performance.
- Adjustable width armrests for comfort and maximized support.
- Lightweight, portable charger for recharging power chair batteries anywhere.
- Optional recline seat upgrade offers a higher seat back and adjustable recline providing additional comfort and support.*
- Optional head rest, oxygen tank, cane, crutches holder, elevated leg rests and many more accessories are available to customize your power chair.*

SPECIFICATIONS

Max weight capacity 300 lbs
 Max speed¹ 4.2 mph
 Max range per charge¹ 12 miles
 Turning radius 24"
 Ground clearance 3.6"
 Maximum step climb 2"

Length 41.9"
 Width 25"

Drive tires (2) 10" x 3.5"
 Caster tires (2) 8" x 2"

Total weight with batteries 176 lbs
 Base only 88 lbs
 Base with batteries 123 lbs

¹ Speed and range will vary with respect to terrain and weight of rider. The Manufacturer reserves the right to alter specifications without notice.

* Optional accessories may not be medically necessary as determined by your doctor or covered by Medicare or third party reimbursement.