




H1000 Dimensions



		SMART SEAT	PERCH**
A	Minimum swivel radius	26.4"	-
	With armrest spacers	27.8"	-
	With extended armrest spacers	28"	-
B	Wall to inner seat back*	3"	3.75"
C	Seat depth	14.25"	12.6"
	With extended seat	16.75"	
D	Minimum folded width to footrest*	15.1"	12"
E	Minimum open width to edge of footrest*	25.6"	27.1"
F	Armrest width - external	21.6"	21.6"
	With armrest spacers	23.25"	
	With extended armrest spacers	25"	
G	Armrest width - internal	18.1"	18.1"
	With armrest spacers	20.25"	
	With extended armrest spacers	22.1"	
H	Seat back height	17.9"	17.9"
I	Footrest to seat height		
	Standard seat option	18.25" - 21.25"	26.5" - 34.5"
J	Minimum footrest height	3.4"	3"
	Minimum track intrusion into staircase	4.5"	4.5"
	Maximum weight limit		
	Standard version	350 lbs	300 lbs
	Extra heavy duty version (XXL)***	440 lbs	

* +1 3/4" for extra heavy duty (XXL) 2-way swivel with extended armrest spacers

** Not available for extra heavy duty (XXL)

*** Not available with hinged rail

Color Choices:

COCOA	SAND	SLATE
RUBY	SAPPHIRE	JADE

Upholstery is made from a wipe clean and flame-retardant material.

All of our stairlifts are tested and approved by worldwide independent inspection services to give you the assurance that our products meet the highest standards in terms of quality and safety.

Ready to take your home back with a Hoveround Stairlift? Call today!

800-625-2190

HOVEROUND
Stairlifts

Hoveround Corporation
2151 Whitfield Industrial Way
Sarasota, FL 34243

800-625-2190
stairlifts@hoveround.com
www.hoveround.com/stairlifts

D80015492

HOVEROUND
Stairlifts



H1000

STRAIGHT STAIRLIFT

H1000

STRAIGHT STAIRLIFT

The H1000 is a heavy duty stairlift that provides minimal intrusion into the staircase. Its anodized aluminum rail offers a stylish and unexposed solution with a hidden gear rack that can safely lift up to 440 lbs. If you need to climb stairs to access your front door or your garden, the H1000's optional outdoor package will give you the freedom to enjoy the outdoors independently again!

- Includes user-friendly retractable seatbelt
Not available on the XXL model
- Each unit comes with 2 handheld remotes to easily retrieve the stairlift
- Prevent unauthorized use with provided key
- Durable add-on packages for commercial and heavy duty applications available
- Limited lifetime warranty on the motor and gear-box; 2 year manufacturer's parts warranty

UP TO
440
POUNDS

Standard weight capacity on the H1000 is 350 lbs., while our H1000XXL boasts the heaviest weight capacity of up to 440 lbs.


A sleek and durable solution that is perfect for both indoors and outdoors

www.hoveround.com/stairlifts



Ergonomic and innovative design features

CONVENIENT



For an easier exit, **swivel** your seat toward the landing once you reach the top, either manually or upgrade to our automatic powered option.

SAFE




Surrounded by **safety sensors**, Handicare stairlifts are designed to stop if meeting an obstruction – keeping you safe and still on your seat.

INTUITIVE




The toggle control was designed in collaboration with occupational therapists to ensure it is **easy-to-use** and operate – simply push with your hand, arm or wrist in the direction you want to go.

VERSATILE



The **outdoor** stairlift is UV protected and weatherproofed to withstand the elements. Includes a waterproof cover so it is always dry and ready.

ADAPTABLE



The **Perch Seat** option is ideal for users with limited hip and knee flexibility. The low profile "perch" style seat means your legs don't have to bend to use it yet you still feel secure.

REVOLUTIONARY



The H1000 can be combined with a powered or manual **hinge** solution to prevent the track from obstructing a doorway at the base of the stairs.

Your first step to independence

For more information call 800-625-2190

HOVEROUND
Stairlifts