




H950 Dimensions



SIMPLICITY SEAT		DIMENSIONS
A	Minimum swivel radius With short armrest	26.75" 25.4"
B	Wall to inner seat back	3"
C	Seat depth	14.5"
D	Minimum folded width to footrest	11.25"
E	Minimum open width to edge of footrest	22.9"
F	Armrest width - external	23.4"
G	Armrest width - internal	18.75"
H	Seat back height	15.75"
I	Footrest to seat height Standard seat option Extended seat post	15.4" - 19.25" 18.1" - 22.1"
J	Minimum footrest height	2"
Minimum track intrusion into stair-case		6.25"
Maximum weight limit Manual hinge or Slide Track		310 lbs 265 lbs



Color Choice: Sand
Upholstery is made from a wipe clean and flame-retardant material.

All of our stairlifts are tested and approved by worldwide independent inspection services to give you the assurance that our products meet the highest standards in terms of quality and safety.

Ready to take your home back with a Hoveround Stairlift? Call today!

800-625-2190



Hoveround Corporation
2151 Whitfield Industrial Way
Sarasota, FL 34243

800-625-2190
stairlifts@hoveround.com
www.hoveround.com/stairlifts

D80015494

HOVEROUND
Stairlifts

H950
STRAIGHT STAIRLIFT



H950

STRAIGHT STAIRLIFT

Our economical, entry level **H950** stairlift is everything but. Not only does the H950 get you safely up and down the stairs, but it is designed to get you there with the highest level of convenience.

- Swivel your seat toward the landing for a safer exit when you reach the top
- Each unit comes with 2 handheld remotes to easily retrieve the stairlift
- Prevent unauthorized use with provided key
- Upholstery provided in a sand-colored wipe clean frame retardent fabric.
- Limited lifetime warranty on the motor and gearbox; 2 year manufacturer's parts warranty

An innovative and economical solution for your straight staircase



www.hoveround.com/stairlifts

Ergonomic and innovative design features

SMART




Stairlifts are battery powered to ensure you are in motion even during a power outage. The **continuous charging** feature charges the battery no matter where it is parked.

SAFE




Surrounded by **safety sensors**, Hoveround stairlifts are designed to stop if meeting an obstruction - keeping you safe and still on

INTUITIVE



The toggle control was designed in collaboration with occupational therapists to ensure it is **easy-to-use** and operate - simply push with your hand, arm or wrist in the direction you want to go.

COMPACT




With a foldable footrest, seat and armrests, the H950 boasts a **slim profile** of only 11.25", ensuring other stair users can still use the staircase comfortably and safely

REVOLUTIONARY



The **Slide Track** moves up the staircase at the same time the seat is moving, leaving the bottom of the stairs completely unobstructed. A manual hinge is also available.

ADAPTABLE



The **Zero Intrusion** option is used when a doorway close to the top of the stairs would be obstructed by the track. *(This option is not available with the hinge track option.)*

Your first step to independence

For more information call 800-625-2190

HOVEROUND
Stairlifts