



The HOVERGlide Lite stair lift by Hoveround is our traditional value leader. It has a quiet drive with soft start and stop which moves the rider effortlessly between levels. The rack and pinion drive system ensures safety and dependability.

specifications

HOVERGlide Lite

- Quiet Drive System with Soft Start and Stop
- DC Drive even runs during power interruptions
- Aesthetic track, finished floor mounts, and integrated charging bases
- Upholstered seat and armrest; folding footrest
- Includes two wireless remote controls
- Pre-mounted on 19" track section
- Service-friendly, two-piece chassis cover design
- Easy access to manual back-up
- 300 lb. weight capacity



Delivered pre-mounted onto a 19" track section for easier, faster installation.

Specifications	
Maximum Weight Capacity	300 lbs. (136.078 kg)
Motor Output Amps at Rated Load	20 amps
Overspeed Governor	No
Track (rail) Type	Extruded Aluminum Natural Finish
Travel	16' Standard; 32' Maximum (9.8 m)
Average Number of Return Trips per Charge	10 - 20
Symmetrical Lift	No
Control in Armrest	Yes
Need to Open Chassis to Install or Change Sides	Yes
Lift Mounts to Steps or Wall	Steps
Minimum Folded Width	15.25" (38.7 cm)
Minimum Footrest Height	8.25" (20.5 cm)
Clear Distance Between Armrests	18" (45.7 cm)
Floor to Seat Height Low / High	23-26" (58 - 66 cm)
Minimum Wall to Stair Side of Rail	11" (27.9 cm)
Seat Depth	14.5" (36.8 cm) / 16" (40.6 cm)
Backrest Height From Top of Seat	13.5" (34.3 cm) / 18" (45.7 cm) (Standard)
Electrical Requirements	100V - 240V AC
Operation Power	24V DC Battery
Speed	19 fpm
Incline Range	25 - 50*
Drive System	Rack & Pinion Gear
Safety Devices	Limit Switches Final Limit Switch Footrest Obstruction Switch Swivel Seat Cut-off Switch Constant Pressure Controls Seat Belt
Warranty:	3-Year on Components 1-Year on Batteries