



PHOENIX HD 3WHL heavy-duty travel scooter



FEATURES & BENEFITS

- Quick and easy disassembly.
- Includes interchangeable panels in Red & Blue.
- Armrests are padded and adjustable.
- Quick connect batteries.
- Ergonomic throttle control.
- Height-adjustable swivel seat with two-tone upholstery.

- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Headlight.
- Swivel seat with fold-down backrest.
- Easy freewheel operation.

SPECIFICATIONS

HCPCS.....	K0801
Weight Capacity	350 lbs.
Top Speed	4 mph
Max Range.....	15 miles
Turning Radius.....	43"
Climbing Angle.....	6°
Ground Clearance.....	3.75"
Dimensions	39" L x 22.5" W
Floor-to-Seat Height.....	20" - 22"
Drive System.....	24V x 350W x 3600rpm
Seat Type.....	Folding
Seat Dimensions	17.5" W x 17" D
Motor.....	24V x 350W x 3600 rpm
Batteries(2)	12V x 20AH

Charger	1.5A Offboard
Brakes.....	Electromagnetic
Front Wheel.....	8" x 2" Flat-free
Rear Wheels	8" x 2" Flat-free
Base Weight63 lbs.
Battery Weight(2)31 lbs.
Seat Weight.....	.21 lbs.
Total Weight	115 lbs.
Heaviest Piece Weight.....	.32 lbs.
Warranty on Frame	Lifetime
Warranty on Electronic	
Controller/Drive Train Components	14 Mths
Warranty on Batteries	6 Mths