



SCOUT 4WHL

compact travel scooter



FEATURES & BENEFITS

- Quick and easy disassembly.
- Interchangeable red and blue panels.
- Armrests are padded and adjustable.
- Quick connect batteries.
- Ergonomic throttle control.
- Height-adjustable swivel seat with fold-down backrest.

- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Delta tiller.
- Easy freewheel operation.

SPECIFICATIONS

Weight Capacity	300 lbs.
Top Speed	4.25 mph
Max Range.....	9 miles
Turning Radius.....	53.75"
Climbing Angle.....	6°
Ground Clearance.....	2.5"
Dimensions	42.5" L x 21.25" W
Floor-to-Seat Height	17" - 19"
Seat Type.....	Folding
Seat Dimensions	16.5" W x 16" D
Controller.....	CTE 45A
Motor.....	24V x 270W x 4700 rpm
Batteries(2)	12V x 12AH

Charger	2A Offboard
Brakes.....	Electromagnetic
Front Wheel.....	8" x 2" Flat-free
Rear Wheels	8" x 2" Flat-free
Base Weight	58 lbs.
Battery Weight(2)	20 lbs.
Seat Weight.....	16 lbs.
Total Weight	94 lbs.
Heaviest Piece Weight.....	39 lbs.
Warranty on Frame	Lifetime
Warranty on Electronic Controller/Drive Train Components	14 Mths
Warranty on Batteries	6 Mths