

# THE *SPITFIRE EX*

COMPACT, EASILY TRANSPORTABLE AND LOADED WITH A NUMBER OF UPGRADED FEATURES

- » Convenient, compact 5 piece design makes tool free assembly and disassembly easy
- » Change the scooter's appearance easily with color change panels in Red, Blue and Silver
- » Two battery options let you select the performance that's best for you
- » 300 lb. weight capacity



SPITFIREXSELL



**drive**<sup>™</sup>

# THE SPITFIRE EX

## FEATURES

- Comfortable, swivel height adjustable seat with stylish two-tone upholstery
- Flip back, adjustable width arms
- 4 mph top speed and a 9 mile maximum cruising range—range can be increased to 15 miles with optional 21AH batteries
- Solid, non-marking, flat-free tires for worry-free travel
- Adjustable angle tiller allows the user to position the controller into the most ideal driving position
- Anti-tip wheels for greater safety
- Seat back folds completely flat, which comes in handy when transporting or storing the scooter
- Includes three sets of interchangeable color panels in Metallic Red, Midnight Blue and Silver
- Comes with easy to remove, quick connect battery pack that can be charged separately from the scooter
- “Finger-Pull” throttle
- Convenient free-wheel release knob
- Armrest-mounted reflectors increase the scooter’s visibility
- Wrap-around delta handles make driving easy, even for those with limited dexterity
- Backlit battery gauge makes it even easier to read
- Standard plastic front basket includes removable liner bag
- Take your Spitfire EX Travel Scooter anywhere you want to go with its simple disassembly feature



## ACCESSORIES

These accessories are sold separately and add functionality and convenience to your power mobility product.

Scooter Arm Tote, Large	S5010
Scooter Arm Tote, Small	S5010SM
Deluxe Scooter Carry Bag	S5000
Cylinder Carry Bag	S5005
Oxygen Cylinder Caddy	SF8010
Scooter Basket	SF8020
Cane Crutch Caddy	SF8030
Sun Shade	SF8050
Oval Bag	AB1100
Carry All Bag	AB1110
Rear Scooter Basket	AB2100
Crutch/Cane Holder	AH1000
Oxygen Cylinder Tank Carrier	AH2000
Drink Holder	AZ0060
Scooter Cover (Small)	AZ1000

## ITEM NUMBERS

### Spitfire EX 1320 (3-Wheel)

SPITFIRE132016FS12	12AH Batteries
SPITFIRE132016FS21	21AH Batteries

### Spitfire EX 1420 (4-Wheel)

SPITFIRE142016FS12	12AH Batteries
SPITFIRE142016FS21	21AH Batteries

## OPTIONS

### Batteries

12AH
21AH

## WARRANTY

Limited Lifetime on frame, 14 months on electronic controller and motor/brake components, six months on batteries.

## SPECIFICATIONS

HCPCS	K0800
Weight Capacity	300 lbs.
Top Speed	5 mph
Maximum Range	9 miles (12AH) 15 miles (21AH)
Turning Radius	46" (1320)   54" (1420)
Climbing Angle	6°
Ground Clearance	3"
Dimensions (L x W)	41" x 20" (1320)   42" x 20" (1420)
Floor to Seat Height	20"-23"
Seat Type	Folding
Seat Dimensions (W x D)	17" x 16"
Controller	PG S-Drive 45A
Motor	24V x 270W x 4700rpm
Batteries (pair)	12V x 12AH   12V x 21AH
Charger	2A Offboard
Brakes	Electromagnetic
Front Wheels	7.5" x 2" Flat Free
Rear Wheels	8.5" x 3" Flat Free
Base Weight	60 lbs. (1320) 66 lbs. (1420)
Battery Weight (pair)	22 lbs. (12AH) 28 lbs. (21AH)
Seat Weight	17 lbs.
Total Weight	1320: 99 lbs. (12AH) 105 lbs. (21AH) 1420: 105 lbs. (12AH) 111 lbs. (21AH)
Heaviest Piece Weight	34 lbs. (1320) 40 lbs. (1420)

SPITFIREXSELL



5 22383 28598 2

**drive**™