




# Smart seat – Handicare Simplicity<sup>+</sup>

## Straight stairlifts



Simplicity <sup>+</sup>	Dimensions
<b>A</b> Minimum swivel radius - With armrest spacers - With extended armrest spacers	25 ¾" 26 ¾" 27 ¾"
<b>B</b> Wall to inner seat back	2 ⅝"
<b>C</b> Seat depth - With extended seat	14 ¼" 16 ¾"
<b>D</b> Minimum folded width to footplate	14 ½"
<b>E</b> Minimum open width to edge footplate	22 ½"
<b>F</b> Armrest width – external - With armrest spacers - With extended armrest spacers	21 ⅝" 23 ¼" 25"
<b>G</b> Armrest width – internal - With armrest spacers - With extended armrest spacers	18 ⅛" 20 ¼" 22 ⅛"
<b>H</b> Seat back height	17 ⅞"
<b>I</b> Footplate to seat height - Standard seat option - Perch seat option	17" - 20 ⅞" 26 ½" - 34 ½"
<b>J</b> Minimum footplate height	2 ¾"
Minimum track intrusion into staircase	7 ¼"
Maximum weight limit - Hinged rail/Sliding Track	310 lbs 265 lbs

Some of the above dimensions depend on the characteristics of the staircase in the home. Please consult the sales manual for further information about these dimensions and features of this seat. All information is provided in good faith but does not constitute a guarantee. Handicare reserves the right to alter designs and specifications without prior notice.

