




Simplicity seat – Handicare 2000 Simplicity

Twin rail curved stairlifts



Simplicity	Dimensions
A Minimum swivel radius - With short armrest	28 3⁄8" 27"
B Wall to inner seat back	4"
C Seat depth	15"
D Minimum folded width to footplate	15 3⁄4"
E Minimum open width to edge footplate	27 1⁄8"
F Armrest width – external	25"
G Armrest width – internal	18 3⁄4"
H Seat back height	16 3⁄8"
I Footplate to seat height*	19 1⁄2" - 22 7⁄8"
J Minimum footplate height	2 3⁄4"
Minimum track intrusion into staircase	4 1⁄2"
Maximum weight limit - Heavy duty version	250 lbs 300 lbs

*Footplate to seat height reduced to 18 3⁄4" if footplate is set to higher position (3 1⁄2").

Some of the above dimensions depend on the characteristics of the staircase in the home. Please consult the sales manual for further information about these dimensions and features of this seat. All information is provided in good faith but does not constitute a guarantee. Handicare reserves the right to alter designs and specifications without prior notice.

